

**MODEL 789A-3 SCAN CONTROLLER**

- Selectable Speeds  
*Local Mode* — 0.1 to 999.9 Å/Min. (Most Monochromators)  
*RS232* — Dependent on Monochromator and Grating (Consult factory.)
- Selectable baud rates: 300, 1200, 2400, 4800, 9600, 19.2K, and 39K (Factory set to 9600)
- Limit interlock system
- Operates in either direction
- Slew function permits high speed go to
- Motor interlock system
- LEDs provide status Information at a glance
- RS232 communications standard, IEEE available
- LED digital display available

The Model 789A-3 scan controller provides the user with a low cost, high resolution, easy to use, stepper motor scanning system.

The Model 789A-3 comes complete with a stepping motor to fit McPherson monochromators or sample handling accessories. The cabling for power and limit switches, controller, and power supply are furnished with this



controller. Selectable features include scanning direction, rate, and motor enable/disable. User friendly front panel controls selection of all functions. Our standard smart serial interface permits selectable features to be changed remotely. No analog to digital or I/O boards required.

For users interested in controlling devices via TTL or C-MOS logic, consult factory for details.

Please specify 110 V/60 Hz or 220/50 Hz with your order, also indicate the model number of the instrument to be fitted if not ordering with the instrument.

**SPECIFICATIONS**

<b>Power</b>	115/220 VAC, 50/60 Hz, 29 Watts
<b>Motor</b>	24 volts @ 1.0 amp/phase (NEMA-23) (Consult factory for other available motors.)
<b>Physical Dimensions</b>	Width: 8.38 In (21.27 cm) Height: 5.22 In (13.26 cm) Depth: 14-38 In (36.51 cm)
<b>Motor Resolution</b>	36,000 step/rev - Standard, Up to 50,000 step/rev available
<b>Weight</b>	8.75 lbs (3.98 Kg)

Command	Description	Command	Description
ESC	ABORT/TERMINATE	S	STORE PARAMETERS AS
@	SOFT STOP		DEFAULTS <sup>2</sup>
^C	RESET*	T	TRIP POINT SET
B	JOG SPEEDS SLOW FAST	+	+ STEP COMMAND
C	RESTORE/INITIALIZE	-	- STEP COMMAND
F	FIND HOME	V	SLEW VELOCITY
G	GO FROM ADDRESS/TRACE <sup>2</sup>	W	WAIT IN MILLISECONDS
I	INITIAL VELOCITY	X	EXAMINE PARAMETERS
J	JUMP TO ADDRESS	Z	DISPLAY POSITION
K	RAMP SLOPE	[	READ NON-VOLATILE MEMORY
L	LOOP ON PORT <sup>2</sup>	\	WRITE TO NON-VOLATILE
M	MOVE AT CONSTANT SPEED <sup>1</sup>		MEMORY
O	ORIGIN SET	]	QUERY HARDWARE STATUS
P	PROGRAM <sup>2</sup>	^	QUERY MOTION
Q	QUERY (list) PROGRAM <sup>2</sup>	_	PROFILE
R	RELATIVE INDEX		

\* all daisy chained devices respond

<sup>1</sup> Continuous Motion (distance not set)

<sup>2</sup> For Internal Program Only (see sample program below)

**Program Example: Utilizing a Model 207 monochromator in conjunction with a 600 G/mm grating.**

**Objective:** To scan the instrument at a rate of 10 Å/Min. over a range of 5 Å. After completion of the scan, change direction and slew back to the starting wavelength at a rate of 2000 Å/Min. Motor resolution is set to 36000 steps/rev.

**Program Entry:**

- P84 (CR)
- 84 V60      60 steps/sec (10 Å/min)
- 87+ 1800    5 Å wavelength range (1800 steps)
- 93 W1      Wait 1 millisecond to change velocity.
- 96 V12000   Change V to 12000 steps/sec. (2000 Å/Min.)
- 98-1800    Change direction & return 1800 steps.
- 102 P (CR)   Stores program, exits to command mode.
  
- Q84 (CR)    To verify or review program.
- (CR)        To scroll through the program.
- G84 (CR)    To run program.